

Glass Fiber Plate DNA Extraction Protocol: Manual Protocol Employing Centrifugation Method

Natalia V. Ivanova, Jeremy R. deWaard & Paul D.N. Hebert


PALL Glass Fiber plate



Bind & wash on a square-well block



Positioning GF plate for elution step
with a PALL collar



July 6, 2006

Glass fiber (GF) plate DNA isolation for recent and archived specimens:

Note: Use **Insect Lysis Buffer** and **PALL2** plate for DNA extraction from Arthropods and centrifuge **PALL2** plates at 6000 g.

1. For 1 plate mix 5 ml of **Vertebrate Lysis Buffer** and 0.5 ml of **Proteinase K**, 20 mg/ml in a sterile container. Add 50 μ l of Lysis Mix to each well of 96-well Eppendorf plate.
2. Add a small amount of tissue (e.g. 2-4 mm of insect leg or 2-3 mm³ of ethanol preserved tissue) to each well of 96-well solid skirted microplate (flame sterilize instruments between samples) cover plate with caps.
3. Incubate at 56 °C for a minimum of 6 hours or overnight to allow digestion.
4. Centrifuge at 1500 g for 15 sec to remove any condensate from the cap strips.
5. Add 100 μ l of **Binding Mix** to each sample using multichannel pipette. Shake vigorously for 10 – 15 sec and centrifuge at 1000 g for 20 sec to remove any sample from the cap strips.
6. Remove cap strips and transfer the lysate (about 150 μ l) from the wells of microplate into the wells of the GF plate (PALL1) placed on top of a square-well block using multichannel pipette. Seal the plate with self-adhering foil.
7. Centrifuge at 5000 g for 5 min to bind DNA to the GF membrane.
8. First wash step: Add 180 μ l of **Protein Wash Buffer (PWB)** to each well of GF plate. Seal with a new cover and centrifuge at 5000 for 2 min.
9. Second wash step: Add 750 μ l of **Wash Buffer (WB)** to each well of the GF plate. Seal with a new self-adhering foil and centrifuge at 5000 for 5 min.
10. For PALL2 plates only (to avoid incomplete Wash Buffer removal): Open the sealing cover, close it and centrifuge the GF plates again for 5 min at 6000 g.
11. Remove the self-adhering foil. Place GF plate on the lid of a tip box. Incubate at 56 °C for 30 min to evaporate residual ethanol.
12. Position a PALL collar on the collection microplate and place the GF plate on top. Dispense 30 – 60 μ l of ddH₂O (prewarmed to 56°C) directly onto the membrane in each well of GF plate and incubate at room temperature for 1 min. Seal plate.
13. Place the assembled plates on a clean square-well block to prevent cracking of the collection plate and centrifuge at 5000 g for 5 min to collect the DNA eluate. Remove the GF plate and discard it.
14. Cover DNA plate with cap strips. DNA can be temporarily stored at 4 °C or at –20 °C for long-term storage.
15. Use 1-5 μ l of the DNA for PCR.

Additional note: square-well blocks could be washed with ELIMINase (or with any other DNA removing detergent), autoclaved and re-used.

Reagents and Stock Solutions for DNA Extraction Using Glass Fiber Plates

Reagents:

Description	Abbreviation	Supplier & Catalogue #
Disodium ethylenediamine tetraacetate · 2H ₂ O	EDTA	Fisher Scientific S311-500
ELIMINase		Decon Labs Inc. 1102
Ethyl alcohol (anhydrous)	EtOH 96%	Commercial Alcohols Inc. 472-06-02
Guanidine thiocyanate	GuSCN	Sigma G9277-500g
Molecular biology grade water	ddH ₂ O	HyClone SH30538.02
Polyethylene glycol sorbitan monolaurate	Tween-20	Fluka 93773
Proteinase K		Promega V3021
Sodium chloride	NaCl	Fisher Scientific S271-3
Sodium dodecyl sulfate	SDS	GibcoBRL 15525-025
Sodium hydroxide	NaOH	Fisher Scientific S318-3
t-Octylphenoxy polyethoxy ethanol	Triton X-100	Sigma T8787-100ML
Tris(hydroxymethyl)aminometane	Trizma base	Sigma T6066-100g
Tris(hydroxymethyl)aminometane hydrochloride	Trizma HCl	Sigma T5941-100g

Disposables & equipment:

Description	Abbreviation	Supplier & Catalogue #
ABGene (Fisher) 8-Strip flat PCR caps	cap strips	Fisher AB-0783
AcroPrep 96 1 ml filter plate with 1.0 µm Glass Fiber media, natural housing	PALL1	PALL 5051
AcroPrep 96 1 ml filter plate with 3.0 µm Glass Fiber media over 0.2 µm Bio-Inert membrane, natural housing	PALL2	PALL 5053
Axyseal sealing film	self-adhering foil	Axygene PCR-SP
Eppendorf twin.tec 96-well microplates	microplate	Fisher E951020427
Matrix Impact2 pipettor, 15 µl-1250 µl, 8-channel	multichannel pipette	Matrix 2004
Matrix Technologies 1250uL Talltip (102mm) Filter tip.		Matrix 8245
PP MASTERBLOCK, 96 Well, 2 ml	square-well block	Greiner 780271
SBS Receiver Plate Collar	PALL collar	PALL 5225

Stock solutions:

Description	Reagents & Weight	Final volume
1M Tris-HCl, pH 8.0	Trizma base	26.5 g
	Trizma HCl	44.4 g
1M Tris-HCl, pH 7.4	Trizma base	9.7 g
	Trizma HCl	66.1 g
0.1M Tris-HCl, pH 6.4	Trizma base	6.06 g
<i>Note:</i> Adjust pH with HCl to 6.4-6.5.		
1M NaCl	NaCl	29.22 g
0.5 M EDTA pH 8.0	EDTA	186.1 g
	NaOH	~20.0 g
<i>Note:</i> Vigorously mix on magnetic stirrer with heater. The disodium salt of EDTA will not go into solution until the pH of the solution is adjusted to ~8.0 by the addition of NaOH. Useful tip – give a brief rinse to NaOH granules with ddH ₂ O in a separate glass before dissolving them.		
Proteinase K (20mg/ml)	Proteinase K	100 mg
<i>Note:</i> Add 5 ml of ddH ₂ O to a 100 mg package of Proteinase K; aliquot by 0.5 ml. Store at – 20°C and do not freeze-thaw aliquots.		

Additional note: Thoroughly wash labware with ELIMINase, rinse with dH₂O. Weigh reagents using a clean spatula, fill up with the molecular grade ddH₂O to the final volume. Filter buffers through 0.2 µm filter into a clean bottle; make smaller volume working aliquots (e.g. 100 ml). Store stock solutions and working aliquots at 4°C.

Working solutions for DNA extraction:

Description	Volume from a stock solution (ml) or weight (g)		Final volume
Vertebrate Lysis Buffer (VLB)			200 ml
100 mM NaCl	1M NaCl	20 ml	
50 mM Tris-HCl, pH 8.0	1M Tris-HCl, pH 8.0	10 ml	
10 mM EDTA, pH 8.0	0.5M EDTA, pH.8.0	4 ml	
0.5% SDS	SDS	1.0 g	
Insect Lysis Buffer			200 ml
700 mM GuSCN	GuSCN	16.5 g	
30 mM EDTA pH 8.0	0.5M EDTA, pH.8.0	12 ml	
30 mM Tris-HCl pH 8.0	1M Tris-HCl, pH 8.0	6 ml	
0.5% Triton X-100	Triton X-100	1 ml	
5% Tween-20	Tween-20	10 ml	
<i>Note:</i> Vigorously mix on magnetic stirrer with heater.			
Binding Buffer (BB)			500 ml
6M GuSCN	GuSCN	354.6 g	
20 mM EDTA pH 8.0	0.5M EDTA, pH.8.0	20 ml	
10 mM Tris-HCl pH 6.4	0.1M Tris-HCl, pH 6.4	50 ml	
4% Triton X-100	Triton X-100	20 ml	
<i>Note:</i> Vigorously mix on magnetic stirrer with heater. If any re-crystallization occurs, pre-warm at 56°C to dissolve before use.			
Wash Buffer (WB)			475 ml
60 % EtOH	EtOH 96%	300 ml	
50 mM NaCl	1M NaCl	23.75 ml	
10 mM Tris-HCl, pH 7.4	1M Tris-HCl, pH 7.4	4.75 ml	
0.5 mM EDTA, pH 8.0	0.5M EDTA, pH 8.0	0.475 ml	
<i>Note:</i> mix well, store at -20°C.			
Binding Mix (BM)			
	Binding Buffer	50 ml	100 ml
	EtOH 96%	50 ml	
<i>Note:</i> stable at room temperature for 1 week.			
Protein Wash Buffer (PWB)			100 ml
	Binding Buffer	26 ml	
	EtOH 96%	70 ml	
<i>Note:</i> stable at room temperature for ~ 1 week, discard if any crystallization occurs.			

Additional note: Weigh the dry components (e.g. SDS or GuSCN) first, then add required volumes of the stock solutions, and fill up with the molecular grade ddH₂O to the final volume. No filtering is required.